

26.9.19

Material za hematologijo W 98051

Platinum White + ABS

$$1) \overset{2}{962} \times \overset{2}{328} \times 18 \quad \overline{\overline{11}} \quad 6 \text{ kom}$$

$$2) \overset{2}{2000} \times \overset{2}{400} \times 18 \quad \overline{\overline{11}} \quad 2 \text{ kom}$$

$$3) \overset{1}{964} \times 390 \times 18 \quad | \quad 7 \text{ kom}$$

$$4) \overset{1}{937} \times 390 \times 18 \quad | \quad 4 \text{ kom}$$

$$5) \overset{1}{964} \times 40 \times 18 \quad | \quad 1 \text{ kom}$$

$$6) \overset{2}{2020} \times \overset{2}{395} \times 18 \quad \overline{\overline{11}} \quad 5 \text{ kom}$$

$$- \text{lesomal} \cdot 1930 \times 980 \times 3 = 2 \text{ kom}$$